

## MBKN3

## NITRILE GLOVE SAFETOUCH® BLACK LITE

Class I Medical Device and Category III Personal Protective Equipment, Powder Free

### FEATURES AND BENEFITS

- Great elasticity.
- Good dexterity.
- Its smooth appearance gives it correct touch and grip of objects.
- Ultra comfortable, this generation allows an excellent stretchability and an optimal comfort.
- Watertight and airtight.

### Recommendations:

- Glove for professional use.
- In the healthcare environment, this glove is used to make barrier against microorganisms between the professional and the patient and their environment.
- In an industrial environment, this glove is intended to protect the wearer from tested chemicals and microorganisms (bacteria, moulds and viruses).
- The powder-free glove reduces the risk of airborne contamination and limits the risk of irritation.

### PRODUCT DETAILS

Product name:	<b>SafeTouch® Black Lite</b>
Product type:	Single use, no sterile
Material:	Nitrile
Colour:	Black
Shape:	Ambidextrous
Cuff:	Rolled
Surface:	Textured fingertips
Powdered:	No
Palm thickness:	0,07 (±0,02mm)
Fingers thickness:	0,10 (± 0,02mm)
Glove length:	≥ 240 mm
Unit weight (M):	3,5g (± 0,2)
Manufacturing country:	Malaysia

Item	Size	
MBKN3-6XS	T 6	XS
MBKN3-7S	T 7	S
MBKN3-8M	T 8	M
MBKN3-9L	T 9	L
MBKN3-10XL	T 10	XL

Marking on packaging:



### INDUSTRIES

Health / Medical / Dental / Industry / Food

### CERTIFICATION AND STANDARDS

Manufactured under an ISO 13485 and ISO 9001 certified quality system.  
 In accordance with requirements of Reg (EU) 2017/745, on medical devices.  
 Comply with the standards EN 455-1/2/3/4 for single use medical gloves.  
 In accordance with requirements of Reg (EU) 2016/425, on personal protective equipments.  
 Comply with the EN 420, EN ISO 374-1/2/4/5, for protective against tested chemicals and microorganisms including viruses.  
 In accordance with requirements of Reg.(EC) 1935/2004 and the French decree of 5th August 2020. This glove is suitable for food contact, with all types of foodstuff.

### PRECAUTIONS BEFORE USE

The glove must be used on intact skin only.  
 The reuse or extended use may lead to infection or cross-contamination.  
 Please refer to your national legislation about disposal of this product once it has been used.

### STORAGE CONDITIONS

Normal conditions of conservation & storage: not to be exposed to moisture and sun, store in a cool environment or at room temperature.  
 Shelf life of the product: 5 years.

### LOGISTIC INFORMATION








#### Specifications outer case



#### Specifications inner case

#### Color code inner and outer

Item	Size mm	Gross Weight kg	QTY/ pallet	Size mm	QTY	Color
MBKN3-6XS		4,1				
MBKN3-7S		4,7				
MBKN3-8M	262 x 252 x 270	5,1	60 cases (5 layers of 12 cases)	245 x 125 x 48	10 dispenser boxes of 100 units	
MBKN3-9L		5,3				
MBKN3-10XL		6,2				

#### EAN Codes

Item	Outer case	Inner case
MBKN3-6XS	3 662 036 020 802	3 662 036 020 857
MBKN3-7S	3 662 036 020 819	3 662 036 020 864
MBKN3-8M	3 662 036 020 826	3 662 036 020 871
MBKN3-9L	3 662 036 020 833	3 662 036 020 888
MBKN3-10XL	3 662 036 020 840	3 662 036 020 895